

TR-ST100BE Wireless Bed Station



TR-ST100BE AURA 2.0 Wireless Bed Station

The TR-ST100BE is a surface mounted bed station with a six foot call cord. A reset switch is located on the front and the call can only be reset at the station. A red LED indicates a call has been placed. Removal of the call cord plug from station also initiates a call. Included are two 3.0V Lithium batteries.

Engineering Specifications:

Model: TR-ST100BE

Dimensions: 5.5" x 4.5" x 1.5"

Power: Two 3.0V Lithium Battery

Battery Replacement: BAT604

Battery Life: 4-5 years